



# Karl Fischer Moisture Titrator

## MKS-510N/MKA-510N/MKS-500



KYOTO ELECTRONICS



## Karl Fischer Moisture Titrator

### MKS-510N and MKA-510N

The Karl Fischer Moisture Titrator MKS-510N and MKA-510N are the result of KEM's many years of experience, which combines the latest technology and advanced engineering with KEM's vast experience in instrumentation, producing finest volumetric Karl Fischer titrator available today.

Karl Fischer titration is the most reliable method for determination of water content. It titrates for quantitative analysis for moisture in solids, liquids and liquid gases. Many of the international standards, such as **ISO, ASTM, DIN, BS, JIS, etc.**, have adopted the Karl Fischer method for moisture determination.

The measurement results are calculated into concentration and necessary data is printed out by the external printer.

For measurement of solid or samples which cannot directly be put into the solvent, the moisture vaporizer ADP-511S works for it. The ADP-511S is easy to operate and maintains steady conditions while vaporizing moisture contained in a sample. The settings of sample boat maneuver, vaporizing temperature and carrier gas running duration, and other conditions for each method are controllable by storing them in memory of the MKS-510N and MKA-510N.

## Features

### Rigid and light resin main casing

Main casing is made of corrosion resistant resin as well as rigid and light in weight.

### Large LCD screen

Large LCD screen shows easy-to-operate dialog messages as well as easy-to-read measurement results including water content and concentration.

### Easy dispensing by pump/stirrer

Built-in dispensing unit is now standard, which pumps solvent in and out of the titration vessel.

### Separated titration vessel

Titration vessel is now separated from main unit so that KF reagent is isolated.

### KF reagent in burette unit

Factor change of KF reagent is now made easy by simplified burette unit.

### Easy operation

The titration condition, calculation formula, printout format, etc. are stored in memory for each individual method so that every method selection assures optimal titration parameters. The printout items can be selected according to user's need.

### Automatic control of moisture evaporator

On-line control of the ADP-511S Evaporator allows solid samples to be analyzed using the same parameters and conditions. The vaporizing temperature, carrier gas running duration, etc. can be stored in memory of MKS-510N and MKA-510N for automatic control.

### Dual-mode titration

MKA-510N consists of dual 10mL direct drive burette as standard. Each burette can work for normal and back titration. Titration using two different factors of Karl Fischer reagents is possible as well. The direct drive burette system will refill by itself if an analysis requires more than 10mL of titrant to find an endpoint.

### Self-diagnostics


The built-in self diagnostic message helps to locate an error or trouble in operation and find solution.

### Interface for external control

The interface for external control to connect Balance, Printer and Personal computer is now standard.

## Principle of Analysis


In Karl Fischer reaction, moisture in sample reacts with Karl Fischer reagent quantitatively. MKS-510N, MKA-510N and MKS-500 make the analysis based on the following formula:



Base: amine, pyridine, etc.

ROH (solvent): 2-methoxyethanol, methanol, etc.

Add extracting solvent to the titration. Titrate moisture from the solvent with Karl Fischer reagent until solvent equilibrium is reached. Add a known amount of sample. Titrate with Karl Fischer reagent having a known factor (mg H<sub>2</sub>O/mL) until the endpoint is reached. The Karl Fischer reagent factor is determined using water in standard or methanol standard. The moisture concentration of the unknown sample is then calculated.



## Applications

The Karl Fischer Moisture Titrators – MKA-510N, MKS-510N and MKS-500 are used for moisture analysis with a variety of natural products, raw materials and industrial products.

### Organic compounds and raw materials:

Organic acid / Alcohol / Ester / Acetal / Ether / Hydrocarbon / Acid anhydride / Acyl chloride / Acid chloride / Nitrogen compound / Halogen compound / Sulphur compound / Peroxide / Carbonyl compound / Hydrate organic salt / Organic acid, etc.

### Inorganic compounds and raw materials:

Hydrate inorganic salt / Inorganic salt / Acid anhydride / Base anhydride / Inorganic acid / Inorganic peroxide, etc.


### Natural products and industrial products:

Medicines / Body tissues / Alkaloid / Capsules / Fertilizer / Agricultural chemicals / Wood / Pulp fibers / Wools / Textiles / Leathers / Cellophane tapes / Synthetic detergents / Soaps / Cosmetics milk / Butter / Cheese / Oils / Fats / Fatty acid / Dehydrated foods / Starches / Grains / Sugars / Caramels / Chocolates / Teas / Coffees / Citric powder / Spices / Gelatin / Seasonings / Alginic acid / Fish meals / Coals / Coal tars / Heavy oils / Petrol / Kerosene / Transformer oils / Lubricants /

Greases / Silicon oils / Fluxes / Benzine / Gases / Liquefied petroleum gases / Freon gases / Vinyl-chloride monomer / Plastics powder / Plastics chip / Ion-exchange resin / Rubbers / Adhesive pigments / Paints / Inks / Dyes / Carbon blacks / Toners / Liquid crystal materials / Photo materials / Ferrites / Metal powders / Desiccants / Ores / Clays / Cement, etc.

### The ASTM standards below have adopted the Karl Fischer method for moisture determination:

ASTM D 1533-96	Standard Test Method for Water in Insulating Liquids (Karl Fischer Reaction Method)
ASTM D 1744-92	Standard Test Method for Water in Liquid Petroleum Products by Karl Fischer Reagent
ASTM D 3277-95	Standard Test Method for Moisture Content of Oil-Impregnated Cellulosic Insulation
ASTM D 4377-93a	Standard Test Method for Water in Crude Oils by Potentiometric Karl Fischer Titration
ASTM E 203-96	Standard Test Method for Water Using Karl Fischer Reagent



**Titration parameter setup**

**Measurement**

**Automatic statistics**

```

TITR.PARA
Method      1
Titr.Mode   1
End Time    30 sec
F.Vol       0.01 ml
T.Speed     Med
Detect.Mode 1
D.Time      0 sec
L.Time      0 sec
I.Time      0 sec
Blank       on
Start       manu
Max.Vol     40 ml
Oven        off
    
```

```

Model :MKA-510N
S/N   :LGB06C13
Sample#
-----
Reagent#
-----
Name#KENTARD
-----
*** R e s u l t ***
Sample No.    01-01
Date 96/12/14 13:18
Wt1          5.6874 g
Wt2          4.7882 g
Net           0.9792 g
Result        0.6418 %
Bur. No.1    3.145 mL
              6.2843 mg
    
```

```

<RESULT>
No. mgH2O  Conc[ mg]
01 10.090   10.090
02 10.086   10.086
03 10.091   10.091
04 10.078   10.078
05 10.082   10.082
06 10.083   10.083
07 10.092   10.092
08 10.102   10.102
09 10.088   10.088
10 10.081   10.081

Statistics
Unit      10
Means    10.086 mg
SD        0.0073 mg
CV        0.0724 %
    
```

## Specification

Type name	MKS-510N	MKA-510N			
Model name	Karl Fischer Moisture Titrator		Piston burette	1) Burette cylinder with piston..... 1 ..... 2	
Measurement method	Volumetric titration method			2) Automatic switching valve for suction/dose..... 1 ..... 2	
Measuring range	1) 0.1mg to 500mg H <sub>2</sub> O 2) 10ppm to 100% H <sub>2</sub> O			3) Backlash mechanism 4) Dispense speed: 0.5mL/sec. 5) Suction speed (2): 20sec/10mL or 80sec/10mL 6) Burette capacity: 10mL, repeat suction/dose to reach preset volume	
Detection range	1) Water content: 0.005 to 500.00mg H <sub>2</sub> O 2) Titration volume: 0.005 to 99.995mL				
Control method	Titration control, EP detection and concentration calculation by microcomputer			Solvent	1) Minimum 30mL (for S-type vessel) 2) Maximum 100mL (for S-type vessel)
Endpoint detection	Two-pin platinum electrode detects liquid resistance compensated polarized potential			Additional function	1) Controls evaporation when ADP-511S is connected 2) Programmable stored settings including heating temperature, aging and measurement sequence
Endpoint wait time	1) 1 to 99 sec. programmable 2) Potential stat which maintains EP status can be set up.		Pump/stirrer	Built-in dispensing mechanism which pumps in/out solvent in titration vessel	
Titration form	Normal titration	1) Normal titration 2) Back titration	External control	RS-232C 1) for Printer 2) for Electronic balance 3) for External computer	
Special functions	1) Titration speed controlled by six steps 2) Automatic start by injecting sample 3) Automatic drift off by dehydrating titration vessel 4) Start time delayed 0 to 9999 sec. 5) Cut off time 0 to 9999 sec. selected 6) Limit time 0 to 9999 sec. to finish titration		Ambient condition	Temperature : 5 to 35 Humidity : less than 85%RH	
Display	1) 240 × 64 dots, 30 digits × 7lines LCD with backlight 2) Displays: (1) Measured water content (2) Processed data (3) Dialog messages: Pre-titration: "Pre-titr" Standby for measurement: "Ready" Stabilized drift: "stable"		Power	AC100 to 120V/200 to 240V, 50/60Hz, 35W	
			Dimension	1) Main unit 274(W) × 458(D) × 575(H)mm 2) Stirrer 118(W) × 225(D) × 320(H)mm 3) Solvent dispenser 240(W) × 170(D) × 280(H)mm	
			Weight	Approx. 12.5kg      Approx. 13.5kg	
Individual method filing	Parameters for normal titration, evaporation by ADP-511S, manual factor measurement by standard, etc. can be stored in five different methods.		Standard components and parts	(1) MKS-510N Main unit.....1      (1) MKA-510N Main unit.....1 (2) Operation manual..... 1      (2) Operation manual..... 1 (3) RS232C Instructions 1      (3) RS232C Instructions 1 (4) Power cord..... 1      (4) Power cord..... 1 (5) Stirrer rotor..... 1      (5) Stirrer rotor..... 1 (6) Stirrer unit..... 1      (6) Stirrer unit..... 1 (7) Stirrer cable..... 1      (7) Stirrer cable..... 1 (8) Tube connector KF 1      (8) Tube connector KF 2 (9) Anti-diffusion nozzle 1      (9) Anti-diffusion nozzle 2	
Indication of endpoint	Electronic beep			Standard accessories	(1) Grease ..... 1 (2) Wrench (8mm)..... 1 (3) Piston rod to pull out head..... 1 (4) Ground wire .....1
Printer	Optional (Recommended printer: IDP-100)				
Calculation	1) Concentration, statistics (mean value, RSD, SD) 2) Recalculation 3) Factor calculation				
Error messages	Erroneous key entry or settings on max. volume, abnormal polar potential or liquid resistance, etc.				

## Optional Accessories

### Evaporator

### ADP-511S

The model ADP-511S Evaporator is used in conjunction with the Karl Fischer Titrator to measure the moisture concentration of plastic pellets or solid samples which are insoluble in Karl Fischer reagents or which contain interfering substance. ADP-511S heats the sample in a closed heating chamber. The vaporized moisture is carried into the titration vessel by nitrogen gas.

#### Features

The magnetic bar moves by remote control the sample from the sample chamber into the oven which eliminates contamination from atmospheric moisture.

A transparent heatproof glass tube permits the sample condition monitored during vaporizing process.

The built-in microprocessor which closely checks the vaporizing condition allows rapid rise and accurate control of heating temperature.

The vaporizing temperature is displayed in three digits for accurate temperature set-up.

The optional external air pump for carrier gas is available for user's convenience.




#### Specification


Type and model name	ADP-511S Moisture Evaporator
Heater	Electrically conductive clear heater glass
Temperature range	Room temp. to 300
Temperature control	1) Control method: Proportional 2) Setting range: 0 to 300 3) Minimum temperature setting: 1 4) Temperature precision: $\pm 2$ 5) Temperature sensor: Chromel-Alumel thermocouple
Display	1) LED digital 2) Temperature display:  3) Flow display:  mL/min
Heater tube	Pyrex glass tube (OD)30 x 270(L)mm
Sample boat	1) Pyrex glass 2) 68(L) x 25(W) x 15(H)mm capacity 16mL
Carrier gas	1) Nitrogen is not included in supplied parts. Nitrogen gas, governor and tubing have to be prepared by user. 2) Air pump is not included in supplied parts.
Gas dryer	1) Silica gel: 100g.....1 2) Zeolite: 100g.....1
Gas flow	100 to 300mL/min
Connection to KF Titrator	1) When connected by Cable #030-3388, MKA-510N/MKS-510N can control ADP-511S 2) For other make of KF titrator, evaporation can be controlled by key entry on ADP-511S
Ambient condition	Temperature : 5 to 35 , Humidity : less than 85%RH
Power	AC100 to 240V, 50/60Hz
Power consumption	150W
Dimension	297(W) x 206(D) x 230(H)mm
Weight	Approx. 7kg
Standard components and parts	(1)ADP-511S Evaporator..... 1 (5) Tube (240mm)..... 2 (2) Heating unit..... 1 (6) Hose joint 6..... 1 (3) Desiccant tube..... 1 (7) Operating manual..... 1 (4) Heater tube..... 1
Standard accessories	(1) Silica gel 500g..... 1 (5) Power fuse 3.15A..... 4 (2) Zeolite 500g..... 1 (6) Cautions sticker..... 1 (3) Bubbler tube..... 1 (7) Sample boat..... 3 (4) Boat pushing rod..... 1

#### ■ Titration vessels

S-type  
(30 ~ 100mL)  
987403001




D-type  
(30 ~ 100mL)  
987403005



N-type  
(50 ~ 150mL)  
987403002



C-type  
(50 ~ 150mL)  
987403003



## Karl Fischer Moisture Titrator

### MKS-500

The MKS-500 is the simple and low cost Karl Fischer Moisture Titrator. Other than ordinary liquid sample, this new titrator can also perform moisture measurement of solid sample that can be dissolved in solvent as well as powder samples. When connected to the optional evaporator, water content of samples like plastic pellet or solid chemicals can also be measured.

#### Features

Moisture measurement by simple and easy key operation.  
By the endpoint detection method of liquid resistance compensation, highly precise measurement of various samples is possible.



Fully GLP/GMP conformed report can be printed out by the optional printer.  
Safety and EMC features conforming to CE marking declaration.

## Specification

Measuring range	1) 0.1mg to 500mg H <sub>2</sub> O , 10ppm to 100% H <sub>2</sub> O 2) 0.005mL to 100mL	Printer	Optional (recommended printer: IDP-100)
Titration form	Normal titration	Calculation	1) Concentration, statistics (mean value, RSD, SD) 2) Recalculation 3) Factor calculation
Display	1) 16 digits x 2 lines LCD with backlight 2) Displays: (1) Measured water content (2) Processed data (3) Dialog messages	External control	RS-232C 1) for Printer 2) for Electronic balance 3) for External computer
Method	4 (Direct, Indirect, Factor, Calib.)	Power	AC100 to 120V/200 to 240V, 50/60Hz, 35W
Solvent	1) Minimum 30mL (for S-type vessel) 2) Maximum 100mL (for S-type vessel)	Dimension	Approx. 280(W) x 450(D) x 480(H)mm
		Weight	Approx. 12.5kg

## Optional Accessories

### Samplers

<p>Bent-type for powder (1 ~ 1.5mL) 987403011</p>	<p>"C" bent-type for powder 987403020</p>	<p>Powder weighing sampler 987403025</p>	<p>Straight-type weighing sampler 987403026</p>	<p>Viscous sampler 987403030</p>	<p>High viscous sampler 987403031</p>
<p>Liquified gas sampler with SUS pressure vessel(100mL) 984333014</p>	<p>Oil sampler 987403013</p>	<p>Spoon-type for Viscous 987403012</p>	<p>Eggplant-type for powder 987403024</p>	<p>Syringe(20mL) 985003225</p>	<p>Syringe inlet port stopper 984333407</p>

**KEM's Karl Fischer Titration range consists of volumetric method ( MKS-510N/MKA-510N/MKS-500 ) and coulometric method ( MKC-510N/MKC-500 ).**

**Coulometric Karl Fischer Titrators :**

**MKC-510N**



**Features of the MKC-510N**

**Titration cell can be increased to 2 units.  
Bromine number and index can be measured.  
Fully conformed to GLP and GMP.  
When the evaporator is connected, vaporizing curve can be displayed.**

Range	10 $\mu$ g - 100mg H <sub>2</sub> O, 89 $\mu$ g - 890mgBr <sub>2</sub>
Sensitivity	0.1 $\mu$ g H <sub>2</sub> O
Display	30 digits X 7 lines LCD with backlight
External I/O	* Printer via RS-232C * Balance via RS-232C * Computer via RS-232C
Output	H <sub>2</sub> O, concentration, dialog messages
Ambient condition	Temperature : 5 - 35 Humidity : less than 85%RH

**MKC-500**



**Features of the MKC-500**

**A low cost titrator, yet offers the same accuracy as the MKC-510N ; 10  $\mu$ gH<sub>2</sub>O measuring range, and 0.1  $\mu$ gH<sub>2</sub>O detection sensitivity.  
Electronic balance, printer and RS-232C external ports as standard.**

Range	10 $\mu$ g - 100mg H <sub>2</sub> O
Sensitivity	0.1 $\mu$ g H <sub>2</sub> O
Display	16 digits X 1 line LCD
External I/O	* Printer via RS-232C * Balance via RS-232C * Computer via RS-232C
Output	H <sub>2</sub> O, concentration, dialog messages
Ambient condition	Temperature : 5 - 35 Humidity : less than 85%RH



**KYOTO ELECTRONICS  
MANUFACTURING CO.,LTD.**

Overseas Division : 8-3 Niban-cho Chiyoda-ku TOKYO 102-0084, JAPAN  
Fax : +81-3-3237-0537, Phone : +81-3-3239-7333

**URL : <http://www.kyoto-kem.com>**

*Distributed by :*



Specifications and design subject to change for improvements without notice. printed in Japan.

311K33C