


This software stores or edits data in memory card when the card which contains measurement results and conditions or sample parameters from KEM's Automatic potentiometric titrator AT-510, Density meter DA-520/DA-510/DA-500 or Refractometer RA-520 is inserted into the PC card drive of your computer. Being Windows[®] 95/98/NT compatible, the resulting data in memory card can be directly downloaded into Microsoft[®] Excel Workbook or stored in CSV format.



(This picture is edited)

■ Feature

- 1 This software can save in personal computer the measurement result, measurement parameter and sample parameter stored in memory card.
- 2 Measurement results from KEM's Automatic potentiometric titrator, Density meter or Refractometer can be exported to Excel[®] workbooks (Excel ver. 7.0/97/2000).
- 3 The output data can also be filed in CSV format. The CSV format data is made best use of in user's Lotus 1-2-3[®] spread sheets.

■ Application

This software can be save in personal computer the method or data files stored in memory card of instrument shown below.

Instrument	File type that can be stored
AT-510	Data file, Method file, Sample file, Calibration record file, Precision check file, APB unit file, Data file (List)
DA-520	Data file, Method file, Sample file, Calibration record file, Precision check file
DA-510	Data file, Method file, Sample file, Calibration record file, Precision check file
DA-500	Data file, Method file, Sample file, Calibration record file, Precision check file
RA-520	Data file, Method file, Sample file, Calibration record file, Precision check file

■ Standard specification

Software	"SOFT-PCR" software in 2 x 3.5 inch Floppy Diskettes
Data storage	1) Export direct to Microsoft Excel 2) Store in CSV format
Filing function	1) All the files in memory card can be saved in computer disk at once. 2) Desired files in memory card can be stored in computer disk. 3) Files in computer disk can be saved in memory card.
Connectable instrument	AT-510, DA-520, DA-510, DA-500, RA-520

■ Required peripherals

Personal computer	Windows 95/98/NT with more than i486DX and Excel 7.0/97/2000
Memory	More than 16MB
Hard disk	Space more than 10MB
Display	More than 640 x 480
Drive	3.5 inch floppy diskette drive and PC card drive
Interface	Printer (Parallel) Port [D-Sub 25 pins]

■ Recommended system

Personal computer	Windows 95/98/NT with more than Pentium 133MHz and Excel 7.0/97/2000
Memory	More than 23MB
Hard disk	Space more than 30MB
Display	More than 1024 x 768
Drive	3.5 inch floppy diskette drive and PC card drive
Interface	Printer (Parallel) Port [D-Sub 25 pins]

KEM
KYOTO ELECTRONICS
 MANUFACTURING CO.,LTD.

Overseas Division : 8-3 Niban-cho Chiyoda-ku TOKYO 102-0084,
 JAPAN Fax : +81-3-3237-0537, Phone : +81-3-3239-7333

URL : <http://www.kyoto-kem.com>

Distributed by :